# Invacare® Flotech® Range

EN	<b>Cushion</b>	<i>User Manual</i>
FR	<b>Coussin</b>	<i>Manuel d'utilisation</i>
DE	<b>Kissen</b>	<i>Gebrauchsanweisung</i>
ES	<b>Cojín</b>	<i>Manual del usuario</i>
IT	<b>Cuscino</b>	<i>Manuale d'uso</i>
NL	<b>Kussen</b>	<i>Gebruiksaanwijzing</i>
PT	<b>Almofada</b>	<i>Manual de instruções</i>
DA	<b>Pude</b>	<i>Brugsanvisning</i>
NO	<b>Pute</b>	<i>Brukerveiledning</i>
SV	<b>Dyna</b>	<i>Bruksanvisning</i>
FI	<b>Tyyny</b>	<i>Käyttöohje</i>

This manual must be given to the user of the product.  
Before using the product, read this manual and save it for future reference!

**Yes, you can.®**

© *Invacare* Ltd

Every effort has been made to ensure that the contents of this publication are fully up-to-date at the time of going to print. As part of its ongoing improvement of products, *Invacare* Ltd reserves the right to modify existing models at any time. *Invacare*® dealers will be notified of any such modifications. Any use of this publication, or of parts thereof, as well as any reproduction of images, must have the written consent of *Invacare* Ltd.

Flotech® is a registered trademark of *Invacare* Ltd.

<b>1</b>	<b>General</b>	<b>4</b>
1.1	General information	4
1.2	Symbols	4
1.3	Warranty	4
1.4	Intended use	5
1.5	Service life	5
1.6	Limitation of Liability	5
<b>2</b>	<b>Safety</b>	<b>6</b>
2.1	Safety information	6
<b>3</b>	<b>Description</b>	<b>6</b>
3.1	Product description	6
<b>4</b>	<b>Usage</b>	<b>7</b>
4.1	Safety information	7
4.2	Using the cushion	7
<b>5</b>	<b>Transport</b>	<b>7</b>
5.1	Safety information	7
<b>6</b>	<b>Maintenance</b>	<b>8</b>
6.1	Inspection	8
6.2	Cleaning and care	8
<b>7</b>	<b>After-use</b>	<b>9</b>
7.1	Storage	9
7.2	Re-use	9
7.3	Disposal	9
<b>8</b>	<b>Technical data</b>	<b>10</b>
8.1	General data	10

# I General

## I.1 General information

Essential nursing care is pivotal in pressure ulcer prevention. This cushion will positively contribute to the outcome of a pressure ulcer prevention care plan.

Education, clinical judgement and action based planning based on vulnerability are fundamental factors in prevention of pressure ulcers.

A range of assessment scales can be used as a formal method of assessing risk from pressure ulcer development, and should be used in conjunction with an informal assessment (informed nursing judgement). Informal assessment is considered to be of greater importance and clinical value.

Please heed all the notes, particularly the safety information, and act accordingly.

Should you require more information on the *Invacare Flotech®* range please refer to the *Invacare Flotech®* Brochure and contact details on the back page of this User Manual.

## I.2 Symbols

### Symbols in this manual

In this User Manual warnings are indicated by symbols. The warning symbols are accompanied by a heading that indicates the severity of the danger.

#### WARNING

Indicates a hazardous situation that could result in serious injury or death if it is not avoided.

#### CAUTION

Indicates a hazardous situation that could result in minor or slight injury if it is not avoided.

#### IMPORTANT

Indicates a hazardous situation that could result in damage to property if it is not avoided.

## Tips and recommendations

Gives useful tips, recommendations and information for efficient, trouble-free use.

## Declaration of conformity

This product complies with Directive 93/42/EEC concerning medical devices. The launch date of this product is stated in the CE declaration of conformity.

## Symbols on the product

	Do not pierce or cut		Do not tumble dry
	Do not put near flame		Machine wash temperature 80 °C (wipedown)
	Do not bleach		Tumble dry Ingress resistant cover only
	Orientation Label		Do not dry clean
	Do not iron		

## I.3 Warranty

We provide a manufacturer's warranty for the product in accordance with our General Terms and Conditions of Business. Warranty claims can only be made through the relevant specialist dealer.

## Standard Invacare Terms

This is to certify that your *Invacare* Flotech® product is warranted by *Invacare* Ltd for a period stated in the Table “Technical Data” of this user guide.

The Warranty of your *Invacare* product is valid from time of shipping. If a defect or fault is discovered the *Invacare* dealer or Local Business Development Manager from whom the appliance was obtained must be notified immediately.

The manufacturer will not accept responsibility for damage caused by misuse or non-observance of the instructions set out in this user guide.

During the period of the warranty any products that have become defective due to faulty workmanship or materials will be renewed without charge.

The warranty will be forfeited should any unauthorized alteration be made to the equipment.

Both warranty and fire retardancy certification will become null and void if non-*Invacare* spares are used on any *Invacare* Flotech® products.

The purchaser's statutory rights under the Consumer Protection Act are not affected.

## Quality and flame retardancy

Quality is fundamental to the company's operation, working within the disciplines of ISO 13485 and ISO9001.

The *Invacare* Flotech® cushions features the CE mark, in compliance with the Medical Device Directive.

The foam and cover used to manufacture the *Invacare* Flotech® range is independently tested and certified in accordance with EN 1021-1 and EN 1021-2.

For further information please contact *Invacare* in your country (addresses see back page of this manual).

## 1.4 Intended use

This pressure reducing cushion is intended to be used as part of an overall pressure ulcer prevention program of care.

This product has been designed to deliver effective pressure reduction and postural support to users, when the product is in normal use which is defined by *Invacare* Ltd as when the Flotech product is placed without additional covers or padding between the user and the support surface.

## 1.5 Service life

We estimate a life expectancy of three years for these products, provided they are used in strict accordance with the intended use as set out in this document and all maintenance and service requirements are met. The estimated life expectancy can be exceeded if the product is carefully used and properly maintained, and provided technical and scientific advances do not result in technical limitations. The life expectancy can also be considerably reduced by extreme or incorrect usage.

The fact that we estimate a life expectancy for these products does not constitute an additional warranty.

## 1.6 Limitation of Liability

The warranty does not extend to the consequential costs from fault clearance, in particular freight and travel costs, loss of earnings, expenses, etc.

- Natural wear and tear.
- Inappropriate or incorrect use.
- Defective assembly or setting-up by the purchaser or third parties.
- Defective or neglectful treatment.
- Use of unsuitable spares.

## 2 Safety

### 2.1 Safety information

#### WARNING

- ▶ Do not use this product or any available optional equipment without first completely reading and understanding these instructions and any additional instructional material such as user manuals, service manuals or instruction sheets supplied with this product or optional equipment. *Invacare* product manuals are available on the internet or your local dealer (Addresses are displayed on the back cover of this manual).
- ▶ If you are unable to understand the warnings, cautions or instructions, please contact a healthcare professional, dealer or technical personnel before attempting to use this equipment – otherwise, injury or damage may occur.

#### WARNING

*Invacare* products are specifically designed and manufactured for use in conjunction with *Invacare* accessories. Accessories designed by other manufacturers have not been tested by *Invacare* and are not recommended for use with *Invacare* products.

The introduction of certain third party products between the cushion surface and the user may reduce or impede the clinical effectiveness of this product.

'Third party products' may include, but are not limited to items including under blankets, plastic sheets and sheepskins, etc.

#### WARNING

Any object between the user and the pressure reducing surface may result in the development of pressure ulcers.

- ▶ Make sure that the support surface in contact with the user is kept free from crumbs and other food debris.
- ▶ Drip cables, stents and other foreign objects must not become entrapped between the user and the pressure reducing surface of the cushion.

#### IMPORTANT

The information contained in this document is subject to change without notice.

- ▶ Check all parts for shipping damage and test before using.
- ▶ In case of damage, do not use.
- ▶ Contact *Invacare*/Carrier for further instructions.

## 3 Description

### 3.1 Product description

The *Invacare* Flotech® range of products consist of cushions, a back support and a lumbar pad.

The cushion range is made of a unique, pressure reducing, ergonomically designed green foam, which is moulded to provide postural support. The foam incorporates a deodorising additive which helps to keep the foam odour free for longer.

The back support has been designed using unique, ergonomically designed green foam, to provide stability to a wheelchair or standard armchair. It can be used in conjunction with any Flotech® cushion.

The Flotech® range also includes the Flotech® lumbar pad, which is designed to be used on chairs or wheelchairs and can improve posture, comfort and support.

The cushion consists of the following parts:

The inner layer consists of foam made from different shapes and densities depending on the model. The Flotech® Plus, Image, Solution, Deep and Solution Xtra cushions also equipped with fluid-Packs for optimization of pressure distribution. The fluid-Packs are composed of medical grade silicone and are non-toxic.

The outer cover is made from an internally coated bi-elastic material and has an anti-slip base with a handle at the front edge of the cushion.

## 4 Usage

### 4.1 Safety information

#### WARNING

- ▶ *Invacare* advise that the patient is measured and assessed by a suitably trained healthcare professional to ensure the correct size and seating requirements to achieve the best pressure and postural care and also position for the individual needs are met.

#### WARNING

It is very important for the patient to reposition themselves, or to be repositioned, on a regular basis. This relieves pressure which helps prevent both tissue compression and potential ulcer formation. This must be based on the clinical judgement of a qualified healthcare professional.

- ▶ Always consult a qualified healthcare professional before using the product.
- ▶ Monitor the patient frequently.

#### CAUTION

- ▶ Please ensure that the cushion is put in place in co-ordination with the orientation label which is situated on the left hand side of the cushion cover at the end of the zip.

#### IMPORTANT

In domestic settings common causes of damage include cigarette burns and the claws of pets that puncture covers, allowing fluid ingress and staining.

#### IMPORTANT

- ▶ To prevent accidental cover damage, please ensure that sharp objects are not placed on the cover. When using patient transfer aids, care should be taken not to damage the Flotech® cover. All transfer aids should be checked for sharp edges or burrs before use.
- ▶ It is important to ensure that the Flotech product is not jammed or damaged by sharp edges when used on wheelchairs.

### 4.2 Using the cushion

1. Remove all packaging before use.
2. Place the *Invacare Flotech*® product directly on the wheelchair or standard armchair.

## 5 Transport

### 5.1 Safety information

#### IMPORTANT

- ▶ Take care when handling the Flotech® product to ensure no damage to the product.
- ▶ Avoid contact with jewellery, nails, abrasive surfaces etc.
- ▶ Avoid contact with walls, door frames, door catches or locks etc.
- ▶ Do not transport in roll cages unless completely protected from the sharp edges of the cage.

## 6 Maintenance

### 6.1 Inspection

Check cushion (foam and cover) for strike-through (this may include fluid ingress, stains, rips or damage) after the release of each patient or on a minimum monthly basis (depending on which occurs first) by a suitably qualified and competent person.

#### Check cushions

1. Unzip the cover completely.
2. Check for any staining on the interior foam.
3. Replace any stained items and dispose of as per local authority procedure.

### 6.2 Cleaning and care

#### ! IMPORTANT

All cleaning agents and disinfectants used must be effective, compatible with one another and must protect the materials they are used to clean. For further information on decontamination in Healthcare Environments, please refer to 'The National Institute for Clinical Excellence' guidelines on Infection Control [www.nice.org.uk/CG002](http://www.nice.org.uk/CG002) and your local infection control policy.

#### Cleaning covers

(Removal of contaminants such as dust and organic matter)

1. Remove cover for laundering.
2. Launder the cover at the recommended temperature between 60 and 80 °C using a diluted detergent solution (instructions on label).

#### ! IMPORTANT

Washing at higher temperatures may cause shrinkage.

#### Drying covers

1. Hang cover from a line or bar and drip dry in a clean environment.

#### ! IMPORTANT

- Tumble drying only for the ingress residence cover as this ensures that the fluorocarbon finish remains effective.
- Tumble dry setting must not exceed 80 °C.
- Do not tumble dry for longer than 10 minutes.
- Dry covers thoroughly before re-fitting to the Flotech product.

#### Disinfecting covers

(Reducing the number of microorganisms)

Please contact your hygiene specialist in the event of contamination.

#### ! IMPORTANT

Ensure that any residual detergent has been removed prior to disinfection.

#### Light soilage

1. Wipe down the cover with a 0.1% Chlorine Solution (1,000 ppm).
2. Rinse the cover with clean water using a single use nonabrasive cloth.
3. Dry the cover thoroughly.

#### Heavy soilage

Where the cover is badly soiled, we recommend cleaning with a dilute cleaning detergent at 80 °C in the washing machine.

Large spillages of blood should be absorbed and removed with paper towels first.

1. Clean up all spillages of bodily fluids i.e. blood, urine, faeces, sputum, wound exudate and all other bodily secretions as soon as possible using a 1% Chlorine Solution (10,000 ppm).
2. Rinse the cover with clean water using a single use nonabrasive cloth, and dry thoroughly.
3. Dry the cover thoroughly.

**! IMPORTANT**

1% Chlorine Solution used on a regular basis can diminish the life of the cover if not rinsed and dried properly.

- Do not use granules.

**WARNING**

- Remove contaminated foams from use.

**CAUTION**

- Keep clear of open heat sources.

**! IMPORTANT**

- Do not use phenols, alcohols, bleaches, or other abrasive materials.

**Replacing covers**

1. Unzip the cover and remove it carefully from the foam core.
2. Place new cover onto the foam core.
3. Then close the zipper.

**! IMPORTANT**

- Ensure that the orientation label on the cover of the cushion matches the orientation printing on the foam core.
- Ensure that the corners are positioned correctly into the corners of the cover.

**Gel spillage**

In the event of an accidental gel spillage from a gel sac, the best way to clean the gel up is:

**On floors**

1. Wipe up the excess gel with a paper towel.
2. Cover the area with a good quality washing up liquid and rub into the area leaving for 30 minutes, to emulsify.
3. Clean area with hot water and a mop then rinse.

4. If required repeat process until spillage has been completely removed.

**On clothes/material**

1. Wipe up excess gel with a paper towel.
2. Rub a good quality washing up liquid into the area and leave for 30 mins.
3. Rinse off using warm water.
4. Wash in washing machine as recommended on the item.

**7 After-use****7.1 Storage****! IMPORTANT**

- Store cushions in a dry environment.
- Store cushions within a protective cover/box.
- Store cushions in a flat state on clean, dry, off-flooring free from sharp edges to avoid any possible damage.
- Never store other items on top of a cushion.
- Do not store cushions next to radiators or other heating devices.

**7.2 Re-use**

The product is suitable for repeated use. The number of times it can be used depends on how often and in which way the product is used.

- Before reuse, clean the product thoroughly, see chapter 6.2 Cleaning and care, page 8.

**7.3 Disposal**

The disposal and recycling of used devices and packaging must comply with the applicable legal regulation.

- Ensure that the cushion is cleaned prior to disposal to avoid any risk of contamination.

## 8 Technical data

### 8.1 General data

Product	Warranty <sup>1)</sup>	Firetesting	Grade ref & Colour	Nominal density range [kg/m <sup>3</sup> ]	Nominal hardness range [N]	Maximum user weight [kg]	Weight of product [kg] <sup>2)</sup>
Flo-tech Lite	3 years	EN 1021-1, 2	FG70/180 green	65 – 70	170 – 190	127	1.5
Flo-tech Lite Visco	3 years	EN 1021-1, 2	LM85/115 blue LM85/150 blue LM85/190 blue	80 – 90	100 – 130 140 – 160 180 – 200	(Soft) 32 (Medium) 51 (Firm) 100	1.26
Flo-tech Contour	3 years	EN 1021-1, 2	FG70/180 green	65 – 70	170 – 190	152 (505 x 480 mm and above) 191	1.38
Flo-tech Contour Lo-Back	3 years	EN 1021-1, 2	FG70/180 green	65 – 70	170 – 190	152	1.38
Flo-tech Contour Visco	3 years	EN 1021-1, 2	LM85/115 blue LM85/150 blue LM85/190 blue LM85/250 blue	80 – 90	100 – 130 140 – 160 180 – 200 240 – 260	(Soft) 51 (Medium) 76 (Firm) 108 (Very firm) 127	1.6
Flo-tech Plus	3 years	EN 1021-1, 2	FG70/180 green	65 – 70	170 – 190	152 (505 x 480 mm and above) 191	2.5
Flo-tech Image	3 years	EN 1021-1, 2	FG70/180 green	65 – 70	170 – 190	140	1.82
Flo-tech Solution	3 years	EN 1021-1, 2	FG70/180 green	65 – 70	170 – 190	140	3.08
Flo-tech Deep Solution	3 years	EN 1021-1, 2	FG70/180 green	65 – 70	170 – 190	140	3.08
Flo-tech Solution Xtra	3 years	EN 1021-1, 2	FG70/180 green	65 – 70	170 – 190	140	3.08
Flo-tech Back Support	3 years	EN 1021-1, 2	FG70/180 green	65 – 70	170 – 190	165	1
Flo-tech Back Support Lumbar	3 years	EN 1021-1, 2	FG70/180 green	65 – 70	170 – 190 80 – 90	165	1.1

Product	Warranty <sup>1)</sup>	Firetesting	Grade ref & Colour	Nominal density range [kg/m <sup>3</sup> ]	Nominal hardness range [N]	Maximum user weight [kg]	Weight of product [kg] <sup>2)</sup>
Flo-tech Lumbar Pad	1 year	EN 1021-1, 2	RX39/200 VC50 055	39 – 41	180 – 200	n/a	n/a

<sup>1)</sup> Warranty is against manufacturing defects.

<sup>2)</sup> Weights can vary depending on size of cushion ordered, average weights used as an indication.

Local requirements / certification available upon request from *Invacare*.

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Invacare A/S, Sdr. Ringvej 37,
Invacare AB, Fagerstagatan 9,
Invacare UK Operations Limited